


## The Spectrum of Hydrogen

### Objective

- To study the emission spectra of the Hydrogen gas.
- Visualize the spectral lines of the Balmer series using a spectrometer.
- Determine Plank's constant,  $h$ .

### Theory

If an electron at a certain energy level in an atom gains sufficient amount of energy, it may jump to higher energy levels and the atom is said to be an excited atom. The atom usually doesn't last for long time in its excited state, but tends to return back to its original state by emitting electromagnetic radiation, which is called photon (see Fig. 1).



**Figure 1:** The emission of photon

If an electron in hydrogen atom, at energy level,  $E_i$ , makes a transition to a lower energy level,  $E_f$ , it will emit a photon with energy,  $E_{ph}$ , equals the difference between the two energy levels, i.e.

$$E_{ph} = E_i - E_f. \quad (1)$$

The energy of the emitted photon is given as

$$E_{ph} = h\nu, \quad (2)$$

where  $h$  is Planck's constant, and  $\nu$  is the frequency. Therefore,

$$h\nu = E_i - E_f. \quad (3)$$

The energy of an energy level,  $n$ , is given as

$$E_n = -\frac{2\pi^2 k^2 m e^4}{h^2 n^2}, \quad (4)$$

and the frequency,  $\nu = c/\lambda$ , where  $c$  is the speed of light, and  $\lambda$  is the wavelength of the photon. Therefore, substituting for  $\nu$ , and  $E_n$  in Eq.(3) gives

$$h^3 = \frac{2\pi^2 k^2 m e^4}{c} \left( \frac{1}{n_f^2} - \frac{1}{n_i^2} \right) \cdot \lambda, \quad (5)$$

where:

$k$  is the electrostatic constant =  $9 \times 10^9 \text{ Nm}^2 \text{ c}^{-2}$ ,

$m$  is the electronic mass =  $9.11 \times 10^{-31} \text{ kg}$ ,

$e$  is the electronic charge =  $1.6 \times 10^{-19} \text{ C}$ ,

$c$  is the speed of light =  $3 \times 10^8 \text{ m/s}$ ,

$h$  is the Planck's constant =  $6.6256 \times 10^{-34} \text{ Js}$ ,

$n_f$  is principle quantum number for the final energy level,

$n_i$  is principle quantum number for the initial energy level.


Therefore, if, for a specific transition,  $n_f$ , and  $n_i$  are known, Planck's constant,  $h$ , can be determined by measuring  $\lambda$  for the emitted photon,.

The wavelength,  $\lambda$ , is measured using a spectrometer (Fig. 2). It consists, mainly, of a diffraction grating, collimator, and a telescope. The collimator directs the light waves emerging through its slit, to the diffraction grating. The telescope, which is free to rotate through various angles, gathers the light waves emerging from the grating and directs them towards a magnifying lens at its end. The diffraction grating disperses (diffracts) the polychromatic light waves at different angles.


The diffraction grating is a piece of glass with equally spaced parallel grooves being ruled on one face. When a parallel beam of polychromatic light falls on one side of the grating, each wavelength is dispersed through a specific angle  $\theta$ , such that

$$n\lambda = d \sin\theta, \quad (6)$$

where  $n$  is order number, and  $d$  is the diffraction spacing constant.



**Figure 2:** The spectrometer



**Figure 3:** Probable transitions in the Hydrogen atom

Figure 3, shows the energy levels of hydrogen atom and the probable transitions for the electron between the levels. For each energy level there is a series of spectral lines.

For instance, all spectral lines produced by electron transitions from any level down to level 1 are called Lyman series. As shown in Fig. 3, there are four series in the hydrogen spectrum. They are Lyman, Balmer, Paschen, and Bracket series. In this experiment, the Balmer series is going to be studied.

## Equipment

- Spectrometer.
- Diffraction grating (300 lines/mm).
- mercury and hydrogen discharge tubes.
- power supply .
- piece of black cloth.

## Procedure

### *Adjustments*

1. Align the telescope with the collimator.
2. Assemble the diffraction grating in the holder.
3. Let the grating be at right angle to the collimator-telescope axis.
4. Place the mercury tube in the discharge socket, and put it in front of the collimator.
5. Cover the area between the collimator and the telescope with the black cloth.
6. Turn on the power to the discharge unit.
7. While looking through telescope, move the light source right and left till you see a clear colored line.
8. Adjust the slit screw to narrow the line seen.
9. Rotate the telescope lens to adjust the plus sign.
10. Try to center the colored lined.
11. Record the value of the scale as your zero value ( $\theta_{\text{zero}}$ ).

### *Determination of d value*

1. Rotate the telescope till you see the first appearance of the green line.
2. Read the angle value,  $\theta'$ , and record in Table I.
3. Calculate the angle  $\theta$ , corresponding to this line, and record in the table.
4. Given the wavelength of the green line of mercury  $\lambda_{\text{HgG}} = 5.46 \times 10^{-7}$  m, determine  $d$ .

- Calculate  $d$  using the value written on the grating, and compare.

### Hydrogen spectrum

- Replace the mercury with the H<sub>2</sub> tube. (use the cloth in holding, while power is off).
- Measure the angles corresponding to the colors of the 1<sup>st</sup> order. Record in Table II.
- Calculate the absolute angle values, and record in the table.
- Calculate  $\lambda$  corresponding to each color line, and record in the table.
- Using  $n_f$  and  $n_i$  values given in the table, Calculate Plank's constant,  $h$ .
- Repeat steps 2-5 for the 2<sup>nd</sup> order. Record in Table III.
- Calculate the average value for  $h$  for Tables II & III, then compare with the typical value.

**Table I** ( $\theta_{zero} = \quad$ )

$\lambda_{Hg}$ (m)	$\theta'$	$\theta$	$d$ (m)	$d_{Theoretical}$ (m)
$5.46 \times 10^{-7}$				

**Table II** (1<sup>st</sup> order colors)

Color	$n_f$	$n_i$	$\theta'$	$\theta$	$\lambda$ (m)	$h$ (J.s)
Faint V.	2	6				
Violet	2	5				
Green	2	4				
Red	2	3				

Avg.  $h =$

**Table III** (2<sup>nd</sup> order colors)

Color	$n_f$	$n_i$	$\theta'$	$\theta$	$\lambda$ (m)	$h$ (J.s)
Faint V.	2	6				
Violet	2	5				
Green	2	4				
Red	2	3				

Avg.  $h =$

## Questions

1. How are the Hydrogen gas atoms excited?
2. What is the purpose of the diffraction grating?
3. What determines the maximum number of the order for a specific wavelength?
4. Is it possible to draw the energy level diagram for the second orbit and above, for the Hydrogen, using your results? How?